



LUX LINE CO

WONDERFLOOR®

HEALTHCARE FLOORING & WALL COVERING

MADE
in
INDIA



Pearl



Mystic Green



Smoke Grey



Mint Green



Ruby Blue



Tangerine



Muddy Beige



Steel Grey



Aqua Blue



Thickness: 2.00 mm
(1.50mm/3.00mm against bulk order)
Roll Width: 1.50 mtr
Roll Length: 20 linear mtr

SIGGMA | Healthcare Floorings


Sigma is specialized homogeneous anti-static flooring specially designed for the healthcare/pharma industry keeping in mind the heavy footfalls, trolley movements etc.

The product has antimicrobial properties throughout its thickness which provides a long time protection against fungal & bacterial growth and is ideal for hospital flooring. The product is specially treated with a "Wonderguard" surface treatment which ensures scratch resistance, resistance to staining, easy maintenance & very easy cleaning of shoe scuff marks.


Salient Features

- Homogeneous, hard wearing
- Random marbled pattern in pastel shades
- Smooth surface for easy cleaning
- Anti bacterial and fungicidal
- **PUR** surface protection for ease in maintenance with stain resistance
- Fire retardant & Low smoke
- The flooring can be made seamless by way of thermal welding
- **Application areas:** hospitals (Pre & post care) blood banks, pathology lab, clean rooms, pharmaceuticals (sterile & production area) food industry etc
- "T" rating for best performance

PUR SURFACE PROTECTION



Wonderfloor's Siggma, Trendo & Wallspro range incorporates **WONDERGUARD** state of the art technology which have exhibited 99.9% inhibition in growth of microbial activity under test conditions and also support clean room condition requirements.



F Green



Ash Wood TW 41



Terracotta



Maple White TW 42




Dk Grey



Red Alm TW 43




D Brown



Eucalyptus TW 44



Icy Blue



Royal Ivory



River Blue



S Grey



Daffodils



Thickness: 2.00 mm
Roll Width: 2 mtr
Roll Length: 15 linear mtr

TRENDO | Multi-Purpose Floorings


Trendo is a state-of-the-art multi-application flooring produced keeping in mind India's quest for the latest trend in comfort, sound & vibration absorption, pleasing shades and a very high abrasion resistance. These floorings are produced with Glass Fibre Impregnation to provide high dimensional stability and minimal shrinkage at joints. The "Wonderguard" surface treatment ensures very low maintenance requirements as dust & grime do not stick to the surface.

Trendo flooring is available in attractive wooden designs & non directional chip design and is most suitable for use in the healthcare industry as it has antibacterial and Anti static property as well. They are also suitable for use in large working areas like offices, libraries, educational institutes, stores, cafeteria and many other similar application.

Salient Features

- Non-directional and wood designs
- Anti bacterial and fungicidal
- **PUR** surface protection for ease in maintenance with stain resistance
- Glass Fibre Impregnation to provide high dimensional stability
- Termites and Moisture resistant
- Reduces noise impact
- High strength wear layer
- Easy to install
- **Application areas:** Hospitals, libraries, educational institutes, stores, cafeterias etc
- **Wear Group P**

Multiple layers for multiple benefits



Wonderfloor's **Sigma, Trendo & Wallspro** range incorporates **WONDERGUARD** state of the art technology which have exhibited 99.9% inhibition in growth of microbial activity under test conditions and also support clean room condition requirements.



Steel Blue



Pearl White



Blue Ice



Rusty Gold




Cappuccino




Slate Grey




Corn




Spring Green




Cloudy Grey




Vanilla



New Grey



Beige



P Green



New Blue



Thickness: 2.00 mm
Roll Width: 2 mtr
Roll Length: 20 linear mtr

ADVENTUS

High Performance
Homogeneous
Flooring

Adventus is a high performance product which can be used in multiple application areas with high footfall. It is easy to maintain. Adventus is non-directional and has a single layered design with multi colour chips.

Salient Features

- Highperformanceflooring
- Homogeneous single layered
- Non-directional design
- PUR surface guard for ease in maintenance and resistance against scratches, scruff marks & stains
- Wear group **T** for superior performance
- Anti bacterial / Anti fungal
- **Application areas:** Hospitals, Pharmaceutical Units, Super markets, Schools, Auditoriums, Library, Large public area, Clean rooms etc.



Thickness: 2.00 mm
Roll Width: 2 mtr
Roll Length: 20 linear mtr


ORBIT

Homogeneous
Conductive
Flooring

Orbit is a homogeneous single layer non-directional conductive flooring with a very heavy qualityof staticcharge control

Salient Features

- Homogeneous singlelayer conductive flooring
- Non directional design
- High quality static control
- Stable and permanent conductivity
- Electrical resistance - 5×10^4 ohms to 10^6 ohms
- Anti-bacterial and fungicidal
- **Application areas:** Operation theater, ICU, MRI, CT Scan, Computer Server Room, Control Room, Laboratories, Blood Banks etc
- Wear Group **P**



Rose White



Ivory



Forest Green



Neo Blue



WALLSPRO installed in Bone & Marrow Hospital Mumbai



Thickness: 1.20 mm
Roll Width: 2.0 mtr
Roll Length: 15 linear mtr

WALLSPRO


Wall Coverings

Wallspro is a wall cladding suitable for use on walls in corridors & common areas of hospitals, labs, clinics, pharmaceutical industry & even in schools. These are having anti-static & anti-fungal properties. These wall coverings are flexible in nature and are impact resistant which provide a long service life once installed. "Wonderguard" surface treatment on Wallspro ensures resistance against stains, which ensures very easy cleaning and easy maintenance.

Salient Features

- Impactresistant
- Antistatic
- Anti bacterial and fungicidal
- Available in a range of soothing pastel shades
- **PUR** surface protection for ease in maintenance with stain resistance
- Glass Fibre Impregnation to provide high dimensional stability
- **Application areas:** Hospitals, dispensaries, pharmaceutical industries, laboratories blood banks, stores, restaurants kitchens, schools, others etc

PUR SURFACE PROTECTION + GFT



Wonderfloor's Siggma, Trendo & Wallspro range incorporates **WONDERGUARD** state of the art technology which have exhibited 99.9% inhibition in growth of microbial activity under test conditions and also support clean room condition requirements.



THOOTHUKUDI GOVERNMENT HOSPITAL
CHENNAI



KK SURGICAL
CHENNAI




ORANGE HOSPITAL
HYDERABAD



SHIFA HOSPITAL
TAMIL NADU



RMD HOSPITAL
PUNE



SDMH HOSPITAL
JAIPUR



OLIVE HEALTHCARE
AHMEDABAD



AMARA HOSPITALS
TIRUPATI



GOVT. HOTEL SALEM
CHENNAI



GHC HOSPITAL
MUMBAI



NATIONAL INSTITUTE OF MEDICAL SCIENCE
JAIPUR



KRS HOSPITAL
CHENNAI

WONDERFLOOR HEALTHCARE CLIENTELE

GOVERNMENT SECTOR	LOCATION
NATIONAL SCIENCE CENTRE	DELHI
AIIMS - DELHI	DELHI
SAFDARJUNG HOSPITAL	DELHI
BABU JAGJIVAN RAM MEMORIAL HOSPITAL	DELHI
SHAH SATNAM JI SPECIALITY HOSPITAL	SIRSA
NARORA ATOMIC POWER STATION	NARORA (PUNJAB)
U.P.S.E.B.	LUCKNOW
PGI	LUCKNOW
I.S.R.O.	AHMEDABAD AND BANGALORE
THERMO FISHER SCIENTIFIC	MUMBAI
BPCL TOWNSHIP	MUMBAI
COOPER MUNICIPAL GENERAL HOSPITAL	MUMBAI
BHABHA ATOMIC RESEARCH CENTRE	MUMBAI
KARNATAKA STATE ROAD TRANSPORT CORPORATION	BANGALORE
POWER GRID CORPORATION OF INDIA LTD	NORTH EAST INDIA
AGARTALA GOVERNMENT MEDICAL COLLEGE	AGARTALA, TRIPURA
ESI HOSPITAL	DURGAPUR, KOLKATA
RIMS	RANCHI
AIIMS	PATNA
CLIENTELE OTHERS	LOCATION
ROCKLAND HOSPITAL	DELHI
ROCKLAND HOSPITAL	MANESAR GURGAON
G.B. PANT HOSPITAL	DELHI
SIR GANGA RAM HOSPITAL	DELHI
METRO HOSPITAL	DELHI / NCR
FORTIS HOSPITAL	NOIDA
IRS HOSPITAL	DELHI
CANCER INSTITUTE & RESEARCH CENTER	LUCKNOW
NIMS COLLEGE JAIPUR	JAIPUR
KAMLA NEHRU HOSPITAL	BHOPAL
VEDANT HOSPITAL	INDORE
JP HOSPITAL	HIMACHAL PRADESH
AHMEDABAD CANCER HOSPITAL	AHMEDABAD
RAJASTHAN HOSPITAL AHMEDABAD	AHMEDABAD
SURAT CIVIL HOSPITAL	SURAT
NOVA EXCEL HOSPITAL	KANPUR
KULVANTI HOSPITAL	LUCKNOW
TENDER HOSPITAL	LUCKNOW
ICON HOSPITAL	LUCKNOW
OLIVE HEALTHCARE	AHMEDABAD
APOLLO HOSPITAL	AHMEDABAD
CADILA HEALTHCARE LIMITED (ZYDUS LIFESCIENCES)	AHMEDABAD
TORRENT PHARMACEUTICALS LTD.	AHMEDABAD
SAKAR HEALTHCARE LTD.	AHMEDABAD
HINDUJA HOSPITAL	MUMBAI
JASLOK HOSPITAL & RESEARCH CENTRE	MUMBAI
S.L. RAHEJA HOSPITAL	MUMBAI
GURU NANAK HOSPITAL	MUMBAI
DEENANATH MANGESHKAR HOSPITAL AND RESEARCH CENTER	PUNE
PRIME SURGICALS	PUNE
RMD HOSPITAL (MANIKCHAND GROUP)	PUNE
NOBLE HOSPITAL & RESEARCH CENTRE	PUNE
YASHODA SUPER SPECIALITY HOSPITAL	HYDERABAD
INDO AMERICAN CANCER HOSPITAL & RESEARCH INSTITUTE	HYDERABAD
NICE CHILDRENS HOSPITAL	HYDERABAD
NARAYANA HRUDAYALAYA HOSPITAL	BANGALORE
AMRIITHA HOSPITAL	BANGALORE
BGS GLENEAGLES GLOBAL HOSPITAL	BANGALORE
FORTIS HOSPITAL	BANGALORE
MIOT INTERNATIONAL HOSPITAL	CHENNAI
APOLLO HOSPITAL	CHENNAI
GOVERNMENT HOSPITAL -SALEM	SALEM
SHIFA HOSPITAL TIRUNELVELI	TAMILNADU
KRS HOSPITAL	CHENNAI
RAMACHANDRA HOSPITAL	CHENNAI
BIG TRUMA HOSPITAL	PATNA
PATNA MEDICAL COLLAGE	PATNA
DHANBAD SADAR HOSPITAL	JHARKHAND
SKMCH HOSPITAL	MUZAFFARPUR
TIBET AYURVEDA	PATNA
CENTRAL INSTITUTE OF PSYCHIATRY	RANCHI
VEERAYATAN EYES & EAR HOSPITAL	RAJGIR - BIHAR

PERFORMANCE CHARACTERISTICS

CHARACTERISTICS	TEST METHOD	SIGGMA	TRENDO	ADVENTUS	ORBIT	WALLSPRO
THICKNESS (mm)	EN 428 / ISO 24346	2.00 (1.50 & 3.00*)	2.00	2.00	2.00	1.20
ROLL WIDTH (mtr)	EN 426 / ISO 24341	1.50	2.00	2.00	2.00	2.00
ROLL LENGTH (mtr)	EN 426 / ISO 24341	20	15	20	20	15
TILE SIZE (mm)	EN 427	300x300/305x305/610x610*	NA	NA	NA	NA
TOTAL WEIGHT (gm/sq mtr) SURFACE	EN 430 / ISO 23997	3200 gm/m ²	3200 gm/m ²	2800 gm/m ²	2800 gm/m ²	1400 gm/m ²
PROTECTION	PUR	PUR PROTECTION	PUR PROTECTION	PUR PROTECTION	NA	PUR PROTECTION
SLIP RESISTANCE	DIN 51130/ASTM D2047	R9	R9	R9	R9	NA
HYGIENE / CLEAN ROOM TEST	ISO 14644-1	Inhibits Growth of Microbes, Class A ≤ 0.10 mm 20 mm MANDREL. PASSES	Inhibits Growth of Microbes, Class A ≤ 0.20 mm 20 mm MANDREL. PASSES	Inhibits Growth of Microbes, Class A ≤ 0.10 mm 20 mm MANDREL. PASSES	CLASS A	Inhibits Growth of Microbes, Class A
RESIDUAL INDENTATION (mm)	EN 433 / ISO 23343-1	≤ 0.40%	≤ 0.10%	≤ 0.40%	≤ 0.10 mm	A
FLEXIBILITY (mm)	EN 435 / ISO 24344	SUITABLE / CONTINUOUS	SUITABLE / CONTINUOUS	SUITABLE / CONTINUOUS	20 mm MANDREL. PASSES	≤ 0.10 mm
DIMENSIONAL STABILITY (%)	EN 434 / ISO 23999	CLASS 1	CLASS 1	CLASS 1	≤ 0.35%	20 mm MANDREL. PASSES
BEARING CASTOR CHAIR	EN 425 / ISO 4918	>99% Inhibits Growth	>99% Inhibits Growth	Inhibits Growth	SUITABLE / CONTINUOUS	≤ 0.10%
SURFACE SPREAD OF FLAME ANTI	BS 476 PART 7	> 8 / PASS BflS1	CflS1	BflS1	CLASS 1	NA
FUNGAL/ANTI-BACTERIAL TEST	ISO846 A,C / ISO 22196	34, 43	31,32	34, 43	Inhibits Growth	CLASS 1
REACTION TO FIRE (kw/m ²)	EN ISO 9239-1	WEAR GROUP T	WEAR GROUP P	WEAR GROUP T	BflS1	>99% Inhibits Growth
PERFORMANCE CLASSIFICATION	EN 685 / ISO 10874	Δ LW = 4db	Δ LW = 14db	Δ LW = 4db	34, 41	CflS1
ABRASION RESISTANCE	EN: 660-2	0, EXCELLENT RESISTANCE	0, EXCELLENT RESISTANCE	EXCELLENT RESISTANCE	WEAR GROUP P	32, 33
SOUND ABSORPTION	ISO 140-8/EN ISO 717-2	> BLUE WOOL STANDARD 5	> BLUE WOOL STANDARD 5	> BLUE WOOL STANDARD 6	Δ LW = 4db	WEAR GROUP M
STAIN RESISTANCE	EN 423 / ISO 26987	SUITABLE	SUITABLE	SUITABLE	AVERAGE RESISTANCE	Δ LW = 2db
COLOUR FASTNESS	ISO 105-B02				> BLUE WOOL STANDARD 6	0, EXCELLENT RESISTANCE >
MARITIME USAGES					SUITABLE	BLUE WOOL STANDARD 5
						SUITABLE

* Against Bulk Order

DIOCTYL
PHTHALATE
FREE

ASTM



MAINTENANCE

Do's: • After laying of flooring, the surface should be cleaned daily with damp mop and weekly or fortnightly with soap solution • Prevent point load indentation by covering your furniture legs with rubber stopper caps • To remove stubborn stains, rub the affected area immediately with a soft coir brush or fine steel wool soaked in a warm soap solution • Scuff marks on the floor can be removed with little application of diesel • Use door coir or suitable mats to minimize accumulation of dust on Wonderfloor vinyl flooring • Use water to remove any chemical spill on flooring.

Don'ts: • Do not use any strong detergent, harsh cleaning agent or dry steel wool • Throwing of lighted cigarette butts or match sticks on the floor should be avoided as it could lead to burn marks on the floor • Do not drag any heavy object or furniture on the flooring to avoid scratches on the floor.

For institutional/bulk enquiries, please e-mail : mkt@rmgpoly.com

or call : +91 9015313059 (Available between 10am-6.30pm)

For export enquiries, please e-mail : info@rmgpoly.com or call : +91 9873138885



CONTACT INFORMATION



Phone Number:
+91 7680075590



Email Address:
info@Luxlineco.com



Website:
www.Luxlineco.com



Location:
8-1-523/300, Brindavan Colony
Tolichowki, Hyderabad,
Telangana - 500008

